



Product Datasheet

H2AFZ (Ab-4) antibody (orb416638)


Description	H2AFZ (Ab-4) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IHC, IP, WB
Immunogen	Peptide sequence around site of Lys (4) derived from Human Histone H2A.Z
Target	H2AFZ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0S5
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot
analysis of K562
whole cel...



Immunoprecipitation
analysis of K562
who...