



Product Datasheet

Acetyl-H2AFZ (K7) antibody (orb416639)


Description	Acetyl-H2AFZ (K7) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, ICC, IF, IP, WB
Immunogen	Peptide sequence around site of Acetyl-Lys (7) derived from Human Histone H2A.Z
Target	H2AFZ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0S5
Dilution Range	ICC:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot analysis of Hela whole cel...



Immunoprecipitation analysis of HepG2 wh...