


Product Datasheet

Formyl-HIST1H1C (K159) antibody (orb416653)

Description	Formyl-HIST1H1C (K159) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Peptide sequence around site of Formyl-Lys (159) derived from Human Histone H1.2
Target	HIST1H1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16403
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cells using Formyl-HIST1H1C (K159) antibody (orb416653).