


Product Datasheet

ASMT antibody (orb416695)

Description	ASMT antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Acetylserotonin O-methyltransferase protein (221-316AA)
Target	ASMT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46597
Dilution Range	IF:1:20-200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of 293 cells...