


Product Datasheet

CRYAA antibody (orb416735)

Description	CRYAA antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Alpha-crystallin A chain protein (98-173AA)
Target	CRYAA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02489
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



IF analysis of
HepG2 cells
lysates
using...