

Product Datasheet

DHX15 antibody (Biotin) (orb416750)

Description	DHX15 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Pre-mRNA-splicing factor ATP-dependent RNA helicase DHX15 protein (371-502AA)
Target	DHX15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43143
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.