

Product Datasheet

GPR1 antibody (Biotin) (orb416778)

Description	GPR1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human G-protein coupled receptor 1 protein (308-355AA)
Target	GPR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46091
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.