



Product Datasheet

HEPACAM antibody (orb416787)

Description	HEPACAM antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Hepatocyte cell adhesion molecule protein (265-415AA)
Target	HEPACAM
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14CZ8
Dilution Range	WB:1:500-5000, IF:1:200-500
Expiration Date	12 months from date of receipt.



IF analysis of MCF-7 cells lysates using...



Western blot analysis of MCF-7 whole cel...