


Product Datasheet

HEXB antibody (orb416791)

Description	HEXB antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Beta-hexosaminidase subunit beta protein (408-528AA)
Target	HEXB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07686
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...