

Product Datasheet

HOXA4 antibody (orb416803)

Description	HOXA4 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Homeobox protein Hox-A4 protein (23-173AA)
Target	HOXA4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00056
Dilution Range	IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.