

Product Datasheet

NRDC antibody (FITC) (orb416877)

Description	NRDC antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Nardilysin protein (32-113AA)
Target	NRDC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43847
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.