



Product Datasheet

PDLIM4 antibody (orb416887)

Description	PDLIM4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Recombinant Human PDZ and LIM domain protein 4 protein (81-243AA)
Target	PDLIM4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P50479
Dilution Range	WB:1:500-5000, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunoprecipitation
analysis of K562
whole cel...



Western blot
analysis of K562
whole cel...