

Product Datasheet

PHOX2A antibody (Biotin) (orb416894)

Description	PHOX2A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Paired mesoderm homeobox protein 2A protein (150-264AA)
Target	PHOX2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14813
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.