

Product Datasheet

CCNO antibody (FITC) (orb41691)

Description	CCNO antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Cyclin-O protein (1-257AA)
Target	CCNO
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22674
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.