


Product Datasheet

RHOD antibody (orb416923)

Description	RHOD antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoD protein (1-80AA)
Target	RHOD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O00212
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of MCF-7 cel...