

Product Datasheet

SLC7A2 antibody (HRP) (orb416948)

Description	SLC7A2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Cationic amino acid transporter 2 protein (604-658AA)
Target	SLC7A2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P52569
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.