

Product Datasheet

FTL3 antibody (HRP) (orb416996)

Description	FTL3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Soybean FT-like protein (1-176AA)
Target	FTL3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	E3NYP3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.