

Product Datasheet

KIAA0319 antibody (HRP) (orb417012)

Description	KIAA0319 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Dyslexia-associated protein KIAA0319 protein (153-268AA)
Target	KIAA0319
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5VV43
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.