

Product Datasheet

gata1-a antibody (HRP) (orb417040)

Description	gata1-a antibody (HRP)
Species/Host	Rabbit
Reactivity	Xenopus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Xenopus laevis GATA-binding factor 1-A protein (1-359AA)
Target	gata1-a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23767
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.