


Product Datasheet

PFKFB4 antibody (orb417063)

Description	PFKFB4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human 6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4 protein (266-368AA)
Target	PFKFB4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16877
Dilution Range	IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...