

Product Datasheet

ZAR1 antibody (FITC) (orb417137)

Description	ZAR1 antibody (FITC)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	FITC
Immunogen	Recombinant Arabidopsis thaliana Receptor protein kinase-like protein ZAR1 protein (167-278AA)
Target	ZAR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9ZU46
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.