

Product Datasheet

SCN10A antibody (orb417203)

Description	SCN10A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Sodium channel protein type 10 subunit alpha protein (992-1099AA)
Target	SCN10A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5Y9
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.