


Product Datasheet

TAX1BP1 antibody (orb417223)

Description	TAX1BP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Tax1-binding protein 1 protein (592-713AA)
Target	TAX1BP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86VP1
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Jurkat whole
C...