

Product Datasheet

Wdr37 antibody (orb417271)

Description	Wdr37 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Mouse WD repeat-containing protein 37 protein (1-496AA)
Target	Wdr37
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8CBE3
Dilution Range	IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.