

Product Datasheet

NLRP7 antibody (FITC) (orb417321)

Description	NLRP7 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 7 protein (89-187AA)
Target	NLRP7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WX94
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.