



Product Datasheet

SCUBE2 antibody (orb417415)

Description	SCUBE2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Signal peptide, CUB and EGF-like domain-containing protein 2 protein (534-703AA)
Target	SCUBE2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NQ36
Dilution Range	WB:1:500-5000, IHC-P:1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of HT29 whole cel...