

Product Datasheet

PNPLA8 antibody (HRP) (orb417440)

Description	PNPLA8 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Calcium-independent phospholipase A2-gamma protein (158-345AA)
Target	PNPLA8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NP80
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.