

Product Datasheet

ONECUT1 antibody (orb417475)

Description	ONECUT1 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Hepatocyte nuclear factor 6 protein (192-292AA)
Target	ONECUT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBC0
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.