

## Product Datasheet

### ONECUT1 antibody (HRP) (orb417476)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | ONECUT1 antibody (HRP)                                            |
| <b>Species/Host</b>        | Rabbit                                                            |
| <b>Reactivity</b>          | Human                                                             |
| <b>Conjugation</b>         | HRP                                                               |
| <b>Tested Applications</b> | ELISA                                                             |
| <b>Immunogen</b>           | Recombinant Human Hepatocyte nuclear factor 6 protein (192-292AA) |
| <b>Target</b>              | ONECUT1                                                           |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol    |
| <b>Form/Appearance</b>     | Liquid                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.     |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | IgG                                                               |
| <b>Clonality</b>           | Polyclonal                                                        |
| <b>Uniprot ID</b>          | <b>Q9UBCO</b>                                                     |
| <b>Dilution Range</b>      | ELISA                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |