

Product Datasheet

EBF3 antibody (orb417483)

Description	EBF3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Transcription factor COE3 protein (416-530AA)
Target	EBF3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H4W6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.