

Product Datasheet

HAO1 antibody (FITC) (orb417545)

Description	HAO1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Hydroxyacid oxidase 1 protein (161-273AA)
Target	HAO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UJM8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.