



Product Datasheet

RAD54L2 antibody (orb417575)

Description	RAD54L2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Helicase ARIP4 protein (1122-1296AA)
Target	RAD54L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y4B4
Dilution Range	WB:1:500-5000, IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Western blot analysis of HeLa whole cel...