

Product Datasheet

TRPV2 antibody (HRP) (orb417588)

Description	TRPV2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily V member 2 protein (1-117AA)
Target	TRPV2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5S1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.