


Product Datasheet

ADAMTS13 antibody (orb417619)

Description	ADAMTS13 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 13 protein (738-905AA)
Target	ADAMTS13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q76LX8
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...