



Product Datasheet

CENPJ antibody (orb417667)

Description	CENPJ antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Centromere protein J protein (1082-1248AA)
Target	CENPJ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HC77
Dilution Range	WB:1:500-5000, IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...



Western blot analysis of K562 whole cel...