



Product Datasheet

EFCAB7 antibody (orb41769)


Description	EFCAB7 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human EF-hand calcium-binding domain-containing protein 7 protein (1-158AA)
Target	EFCAB7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A8K855
Dilution Range	WB:1:2000-10000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...



Immunohistochemical staining of human li...



Immunofluorescence analysis of HepG2 cel...