

Product Datasheet

EPM2A antibody (FITC) (orb41783)

Description	EPM2A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Laforin, isoform 9 protein (1-331AA)
Target	EPM2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B3EWF7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.