


Product Datasheet

SFRP5 antibody (orb417927)

Description	SFRP5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Secreted frizzled-related protein 5 protein (169-261AA)
Target	SFRP5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5T4F7
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...