


Product Datasheet

SRGAP3 antibody (orb417959)

Description	SRGAP3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human SLIT-ROBO Rho GTPase-activating protein 3 protein (709-955AA)
Target	SRGAP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43295
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of A549 cell...