

Product Datasheet

TFB1M antibody (FITC) (orb417969)

Description	TFB1M antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Dimethyladenosine transferase 1, mitochondrial protein (59-194AA)
Target	TFB1M
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WVM0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.