

Product Datasheet

MGAT3 antibody (orb418031)

Description	MGAT3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Beta-1,4-mannosyl-glycoprotein 4-beta-N-acetylglucosaminyltransferase protein (435-533AA)
Target	MGAT3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q09327
Dilution Range	IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.