

Product Datasheet

beta-mammal toxin Css4 antibody (orb418059)

Description	beta-mammal toxin Css4 antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Centruroides suffusus suffusus Beta-mammal toxin Css4 protein (1-66AA)
Target	Beta-mammal toxin Css4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60266
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.