

Product Datasheet

RPE65 antibody (orb418143)

Description	RPE65 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Retinoid isomerohydrolase protein (348-484AA)
Target	RPE65
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16518
Dilution Range	IF:1:200-500
Expiration Date	12 months from date of receipt.