

Product Datasheet

PAG2 antibody (orb418147)

Description	PAG2 antibody
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Bovine Pregnancy-associated glycoprotein 2 protein (22-376AA)
Target	PAG2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q28057
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.