



Product Datasheet

NADK2 antibody (orb418151)


Description	NADK2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human NAD kinase 2, mitochondrial protein (16-156AA)
Target	NADK2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q4G0N4
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of HL60 whole cel...



Immunohistochemistry of paraffin-embedde...