



Product Datasheet

PNPLA5 antibody (orb418215)


Description	PNPLA5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Patatin-like phospholipase domain-containing protein 5 protein (90-205AA)
Target	PNPLA5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z6Z6
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Colo320 whole ...



Immunofluorescence staining of HepG2 cel...