

Product Datasheet

PNPLA5 antibody (Biotin) (orb418218)

Description	PNPLA5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Patatin-like phospholipase domain-containing protein 5 protein (90-205AA)
Target	PNPLA5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z6Z6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.