



Product Datasheet

LAPTM4B antibody (orb418219)


Description	LAPTM4B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Lysosomal-associated transmembrane protein 4B protein (174-226AA)
Target	LAPTM4B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86VI4
Dilution Range	WB: 1:500-5000, IHC-P: 1:500-1000, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western Blot analysis of MCF-7 whole cel...



Immunohistochemical staining of human lu...



Immunohistochemical staining of human li...