

Product Datasheet

G6PD antibody (FITC) (orb41823)

Description	G6PD antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Glucose-6-phosphate 1-dehydrogenase protein (80-185AA)
Target	G6PD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11413
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.