

Product Datasheet

LGI4 antibody (HRP) (orb418256)

Description	LGI4 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Leucine-rich repeat LGI family member 4 protein (292-412AA)
Target	LGI4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N135
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.